




Hydromassage unit

REF. 1098.00



DESCRIPTION :

The hydromassage unit includes a control panel placed at the edge of the pool or fixed on a stainless steel support (option 1098.10), an overpressure group placed in the technical room and the piping group.

OPERATION :

The suction of the pump is connected to a Ø50mm valve, needed, on the suction of the filtration group.

The discharge of the pump is carried out through the lance of massage, the pressure being regulated by a control valve.

SECURITY :

Automatic pump shutdown after 30 minutes of continuous operation.

FEATURES :

Polyester panel equipped with the following equipment :

- Pressure gauge dial.
- Pressure variation valve.
- Branch pipe valve.
- Branch pipe holder.
- Branch pipe and its handle nozzle holder(s).

Overpressure group :

- Type : multicellular pump
- Flow : 3,5m³/h
- Pressure : 0 - 10 bars (adjustable)
- Power : 1,1 KW
- Velocity : 2,835 Tr/mn
- Current : 400 V 50 Hz threephase +neutral+ ground

OPTION :

1098.10 Stainless steel support.

TECHNICAL NEEDS :

To see the technical needs sheet.