

The natural, economic and ecological destruction of waste.

Macerator to collect disposable items

REF. 2103.00

Macerator

Reduce the sources of infection, put the caregiver back at the heart of his profession, offer a comfortable and hygienic solution for the patient: the macerator is the answer to disposal of human waste and single use medical pulp containers.





Hygiene



UNIQUE FEATURES







- Undeniable benefit in terms of hygiene.
- · Limit the risk of cross-contamination.
- Simplification of basin rotation protocols with easy management of clean and dirt.
- Optimises the working time of nursing staff.
- Very quiet.
- Environmentally friendly products (HQE standards).
- Economic use: 15 to 22 litres of water per cycle depending on model and without chemicals.
- Wide variety of single use items and accessories.




The macerator consists of pulverisation and disposal unit and a range of single-use items made from natural cellulose. The system, based on a short cycle, optimizes operating costs while fully respecting the environment.

Macerator to collect disposable items

REF. 2103.00

	QUATTRO VANGUARD
Capacity	4 Items
Cycle Time	Configurable between 82 & 180 seconds
Water consumption per cycle (L)	Between 15 & 22 Litres
Power	500W
Consumption per cycle	0.01KWh
Voltage/Phases/Frequency	240V/1ph/50Hz
Method of opening/closing	Hands free foot control
Odour control	Natural siphon
Start-Up function	Automatic
Function/Dialogue display	LED
Unloaded weight	58.5Kg +/-2kg
Access for maintenance	Front and Top
Floor area	524mm x 411mm
Particle size	2mm - 6mm
Grinding Method	Control
Internal surfaces	Antimicrobial
Noise Level	54dB +/-2dB (In night mode)
Typical Lifespan	15 Years
Integrated disinfection device	Yes
* 4 Items = 4 Standard basins or 4	urinals or 2 large basins

Quattro Vanguard



${\bf ACCESSORIES} \ {\bf FOR} \ {\bf MACERATOR} \ {\bf QUATTRO} \ {\bf VANGUARD};$

- 2103.10 Connection kit
- Raisers available (contact us)

